

Manufacture	Model	Manufacture Magnetron Number	Teledyne e2v Magnetron Number	Manufacture Thyatron Number	Teledyne e2v Thyatron Number
ACCURAY 	CyberKnife M6 version			CX1159G	CX1159
	CyberKnife all except M6 version			8503B CX1191D	8503C -
ELEKTA    	Axesse	MG6370(E)	MG6028	CX1154L	-
		MG6364	MG6030		
	Compact	MG6250, MG6165, MG6493, MG5193	MG5193	CX1551	CX1551G
		MG6160-06, MG6205, MG6405	MG6062		
	Infinity	MG6370(E)	MG6028	CX1154L	-
		MG6364	MG6030		
	Unity	MG7095, MG7195, MG7395, MG7495	MG7095	CX1154L	-
		MG6364	MG6030		
	Precise	MG6370(E)	MG6028	CX1154L	-
		MG6364	MG6030		
	SL Series (Fast tuned)	MG6370(E)	MG6028	CX1154L	-
	Sinergy	MG6370(E)	MG6028	CX1154L	-
		MG6364	MG6030		
	Versa HD	MG6370(E)	MG6028	CX1154L	-
		MG6364	MG6030		
	SL48	MG6330	MG5125	CX1140LD, CX1140LE	CX1150LG
SL75-5	MG6330	MG5125	CX1140GA	CX1140A	
SL75-10, SL75-14 & -15			CX1140GA	CX1140A	
Mechanically tuned magnetrons * not for electromagnetically fast tuned systems.	MG6363	MG5028	CX1140GA	CX1140	

Manufacture	Model	Manufacture Magnetron Number	Te2V Magnetron Number	Manufacture Thyratron Number	Te2V Thyratron Number	
<p>VARIAN</p> 	Clinac 600CD	MG5349, MG6190, MG6290, MG6490	MG6090	8503AF, 8503AG	8503K	
		MG6195, MG6253, MG6453	MG6053			
	Clinac 600C	MG6160-06, MG6205, MG6405	MG6062	8503AF, 8503AG	8503K	
					CX1159	CX1159G
					CX1159	CX1159G
	Clinac 600SR	MG6250, MG6165, MG6493, MG5193	MG5193	8503AF, 8503AG	8503K	
					CX1159	CX1159G
		Clinac 6EX	MG6195, MG6253, MG6453	MG6053	503AF, 8503AG	8503K
					CX1159	CX1159G
		Clinac 6-100C			8503AF, 8503AG	8503K
8503AF, 8503AG						8503K
	Clinac 1800C, 2100C, 2300C			8503AF, 8503AG	8503K	
				CX1140LE, CX1140LD	CX1150LG	
				8503AF, 8503AG	8503K	
	Clinac 21EX, 23EX, IX/CX			8503AF, 8503AG	8503K	
				CX1140LE, CX1140LD	CX1150LG	
	Trilogy, Silhouette			8503AF, 8503AG	8503K	
				CX1140LE, CX1140LD	CX1150LG	
	Unique	MG5349, MG6190, MG6290, MG6490	MG6090	CX1159	CX1159G	
		MG6195, G6253, MG6453	MG6053			CX1140LE, CX1140LD
	TrueBeam (& VitalBeam, Edge)				8503AF, 8503AG	8503K
CX1140LE, CX1140LD					CX1150LG	

MATRIZ DE EQUIVALÊNCIA DE PRODUTOS
Magnetrons e Thyratrons

Manufacture	Model	Manufacture Magnetron Number	Te2V Magnetron Number	Manufacture Thyatron Number	Te2V Thyatron Number	
SIEMENS 	Mevatron 6-67	MG6220	MG5125	CX1159	CX1159G	
	Mevatron 6-67, 6000/7000 series, 4-10			CX1159	CX1159G	
	Mevatron 4-10	MG5349, MG6190, MG6290, MG6490	MG6090			
		MG6195, MG6253, MG6453	MG6053			
	MD				CX1159	CX1159G
	ME	MG5349, MG6190, MG6290, MG6490	MG6090		CX1159	CX1159G
		MG6195, MG6253, MG6453	MG6053			
Primus	MG6250, MG6165, MG6493, MG5193	MG5193		CX1159	-	
	MG6160-06, MG6205, MG6405	MG6062				
TOMO 	Radixact	MG5349, MG6190, MG6290, MG6490	MG6090			
		MG6195, MG6253, MG6453	MG6053			
	H Series	MG6250, MG6165, MG6493, MG5193	MG5193			
		MG6160-06, MG6205, MG6405	MG6062			
		MG7095, MG7195, MG7395, MG7495	MG7095			
VIEWRAY	MRIdian			CX1154L	-	
ZAP	ZAP Surgical	MG7095, MG7195, MG7395, MG7495	MG7095			
VAREX	Linatron	MG6160-06, MG6205, MG6405	MG6062			
VERO	C Bend			CX1528	-	
UNITED IMAGING		MG7095, MG7195, MG7395, MG7495	MG7095			
SHINVA	XHA600D, XHA1400					

Magnetrons e Thyratrons

Manufacture	Model	Manufacture Magnetron Number	Te2V Magnetron Number	Manufacture Thyatron Number	Te2V Thyatron Number
SHINVA	XHA600D, XHA1400	MG5193	MG5193	CX1159	CX1159G
		MG6160-06, MG6205, MG6405	MG6062		
		MG6250, MG6165, MG6493, MG5193	MG5193		
REFLEXION	ZAP Surgical	MG7095, MG7195, MG7395, MG7495	MG7095		